


Recombinant *Sulfolobus islandicus* Chromatin protein Cren7 (creN7)

Catalog Number: CSB-EP502491FPB



| | |
|------------------------|---|
| Product Name: | Recombinant <i>Sulfolobus islandicus</i> Chromatin protein Cren7 (creN7) |
| Catalog Number: | CSB-EP502491FPB |
| Relevance : | A probable chromatin protein, binds double-strand DNA without sequence specificity. Constrains negative DNA supercoils. |
| Mol. Weight: | 36kDa |
| Product Info : | GST-tagged |
| Source: | <i>E. coli</i> derived |
| Image: |  |
| Purity: | >90%(SDS-PAGE) |
| Storage Buffer: | PBS, pH 7.4, 50% glycerol |
| Storage : | Store at -20°C, for extended storage, conserve at -20°C or -80°C. |
| Notes : | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| AA sequence: | MSSGKKAVKVKTPAGKEAELVPEKWWALAPKGRKGVKIGLFKDPETGKYFRHKLPDDYPI |
| References: | "Biogeography of the <i>Sulfolobus islandicus</i> pan-genome." Reno M.L., Held N.L., Fields C.J., Burke P.V., Whitaker R.J. Proc. Natl. Acad. Sci. U.S.A. 106:8605-8610(2009) |