




FAM53B Antibody

Product Code	CSB-PA617903LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14153
Immunogen	Recombinant Human Protein FAM53B protein (1-305AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:2000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein FAM53B (Protein simplet), FAM53B, KIAA0140 SMP
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	FAM53B

Image



Western blot

All lanes: FAM53B antibody at 2µg/ml

Lane 1: MCF-7 whole cell lysate

Lane 2: A549 whole cell lysate


Lane 3: HCT116 whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 46, 34 kDa

Observed band size: 46 kDa



Immunohistochemistry of paraffin-embedded

human breast cancer using CSB-

PA617903LA01HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.